BenCorpp[®]
Improving patients' life

**Why people need
*Trach-HME?***


**SpeVaX[®] HME with
built-in Speech Valve**

**HME is indicated in
patients breathing
spontaneously via an
ET tube or a
tracheostomy tube**

When an **ET tube is used?**

**ET intubation is used to provide
airway protection**

**It can be performed during
general anaesthesia for a
surgical procedure or it can be
necessary for a patient with
respiratory insufficiency**

Patients with an ET tube experience the following diseases:


**Normal Heat and Moisture
function is bypassed**

**Exposure to dry gases
causes a thickened mucus**

**Decreased ciliary action
and increased secretions**

Increased risk of infection

Inability to speak

...and here comes

SpeVaX[®]

Small and light weight


- **Minimizes drag on the tracheostomy or tracheal tube**
- **Fits onto any trach hub**
- **Up to 24 hour use**

Incorporated heat and moisture properties

- **Eliminate the need for active heat and humidifying resources**
- **Decrease mucosal drying**
- **Enhance patient comfort and hygiene**
- **Allow for increased patient mobility**

Integrated Suction Valve

- **Facilitates suctioning**
- **Allows inspiration or exhalation**
- **Eliminates the need for device removal when suctioning**
- **Lowers risk of infection**
- **Eliminates the need for a pressure relief valve**


Why **SpeVaX**[®]?

**Digitally Occluded
Speech Valve
Improves speech
and hygiene**

**Latex Free
Meets hospital
requirements**

SpevaX®

Trach-HME with built-in Speech Valve

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BenCorpp®
Improving patients' life

SpevaX®

Trach-HME with built-in Speech Valve

SpeVaX® is a tracheal heat and moisture exchanger with built-in speech valve indicated in patients breathing spontaneously via an ET tube or a tracheostomy tube. It is recommended for patients weighing no less than 10 kg.

Recommended tidal volume range	50 - 1000 ml
Compressible volume	9,5 ml (dead space)
Moisture output	VT= 500 ml: 13.1 mg/l
Resistance to flow	31 ± 5 Pa a 30 l/min
Weight	2,9 g
Length	22 mm
Height	28 mm
Width	38 mm
Non-conductive	

Features

SpeVaX® should be connected to the open end of tracheostomy tube, and optimal performance is achieved after a few breaths.
SpeVaX® should be replaced every 24 hours, or more often if necessary, to prevent resistance due to the accumulation of secretions.
Before use, check the integrity of the product and valve by closing and reopening it.
Store in a clean, dry, dark place at room temperature.
SpeVaX® is equipped with a spring valve which can be easily occluded with a finger to facilitate speech. After removing your finger, the valve will automatically reopen.

Contraindications

Do not use the valve exceeding the recommended tidal volume range, as dead space can cause CO₂ retention in case of too low tidal volumes, or insufficient humidification in case of very high tidal volumes.
Do not use SpeVaX® on dehydrated patients with heavy secretions from the lungs and airways, or in combination with external humidifiers or nebulizers.

Caution

Single use product: do not sterilize or reuse.
Discard after use.

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BenCorpp Srl