
SpeVaX® – Istruzioni per l'uso

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SpeVaX® – Instructions For Use

Trach-HME with Built-in Speech Valve

Manufacturer: **BenCorpp® Srl**
Via Imperiale, 279
IT - 41037 – Mirandola (MO)
info@bencorpp.com

Indications

Heat & moisture exchanger for patients breathing spontaneously via an ET tube or a tracheostomy tube in the hospital or at home. SpeVaX is recommended for patients weighing no less than 10 kg.

Technical specification

Recommended tidal volume range	50 - 1000 ml
Compressible volume	9,5 ml (dead space)
Moisture output	V _T = 500ml:13.1 mg/l
Resistance to flow	31 ± 5 Pa a 30 l/min
Weight	2,9 g
Length	22 mm
Height	28 mm
Width	38 mm
Non-conductive	

Connections

15 mm female connector for ET tube or tracheostomy tube. SpeVaX® has a valve for forced inspiration/expiration or suctioning.

Instructions for use

The device should be placed at the open end of the ET tube or tracheostomy tube. Check that the connection is properly fixed. Normal efficiency is reached after a few breaths. Should be replaced every 24 hours or more frequently, to avoid increase in resistance due to accumulation of secretions.

The maximum continuous use time of the device must be less than 30 days. In the event that continuous use is necessary for a period exceeding 30 days, it is necessary to seek medical advice.

Before use, check product for deformity and valve function by closing the valve and releasing it. Do not use if the device shows any signs of defect. Do not use if the package is opened or damaged.

Store in clean, dry, dark, temperate conditions.

Speech valve

SpeVaX® contains a valve with a spring that can easily be occluded with a finger facilitating speech. After releasing the finger, the valve will open automatically.

Contraindications

Do not use beyond recommended tidal volume range, as the added dead space may cause CO2 retention at too low tidal volumes. A too high tidal volume may lead to unsatisfactory humidification. Do not use on dehydrated patients or patients with very heavy secretions from the lungs and airways.

Caution

Non-sterile, single use only: do not sterilize or reuse, as the intended function will be affected, and this might increase the risk of possible transmission of infection. Do not use the product after the expiration date. The expiry date refers to the intact and correctly stored product.

Destroy after use. Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

In case of undesirable effects, discontinue use, consult your doctor and inform the manufacturer.

In the event of serious accidents occurring related to the use of the device SpeVaX® Trach-HME with Built-in Speech Valve, send a report to the manufacturer (info@bencorpp.com), to the Italian Ministry of Health or to the competent authority of the EU Member State in which the product has been used.

GRAPHIC SYMBOLS CONTAINED IN DEVICE LABELING

Consult Instructions For Use

Cautions on instructions for use

Manufacturer

Medical Device

Batch Code

Do Not Re-use

Catalog Number

Keep Away from Sunlight

Use By

Keep dry

Unique Device Identifier

Do not use if package is damaged