


DilatationX[®]



Dilation Balloon System


BenCorpp[®]
Improving patients' life

DilataX[®]


Dilation Balloon System

Features

- Indicated for upper airway dilation
- Atraumatic access
- Low profile, high pressure and non-compliant balloon for narrow stenosis
- Preloaded stylet for a faster and safer procedure
- Marker bands to ease balloon placement on stenotic area



Balloon placed across stricture



Balloon inflation dilates stricture

CATALOG NUMBER	DIAMETER (mm)	BALLOON LENGTH (mm)	CATHETER LENGTH (mm)	MAX PRESSURE	INFLATION VOLUME
TLA 005-30-500	5	30	500	16 atm	6-8 cc
TLA 007-30-500	7	30	500	16 atm	6-8 cc
TLA 010-40-500	10	40	500	12 atm	12 cc
TLA 012-40-500	12	40	500	10 atm	12 cc
TLA 014-40-500	14	40	500	10 atm	12 cc
TLA 016-40-500	16	40	500	8 atm	12 cc
TLA 018-40-450	18	40	450	8 atm	12 cc